



入門指南 - USB-to-HDMI 轉換器 DFG/USBtoHDMI

典型的機器視覺應用使用專用的電腦進行影像處理和分析。然而，某些簡單的監控和可視化任務不需要專屬級別的影像處理。這類任務最有效的解決方法是直接地從相機傳送影像到顯示器，從而免除了對專用電腦硬體的需求和相關成本。USB-to-HDMI轉換器即專為此目的而設計，可從The Imaging Source的USB3相機影像數據生成HDMI信號。諸如曝光和增益等相機設置可以透過螢幕選單進行配置，並於啟動時自動恢復。

入門指南

請執行以下步驟進行轉換器設置：

- 將HDMI顯示器連接至轉換器的微型HDMI端口
- 連接電源插頭(USB-C)
- 若尚未連接至相機，請相連接至USB3端口(藍色)
- 若需更改相機設置，請將標準USB鍵盤連接到設備，然後按Enter打開主選單

產品特色

- 支援所有The Imaging Source DMK/DFK 33, 37 及 38 系列相機
- 解析度高達4K(3840x2160)，其幀速率達15fps；全高清時，幀速率高達60fps
- 支援顯示器高達4K原始解析度
- 易於使用的螢幕選單，可透過USB鍵盤進行相機和顯示設置

供貨範疇

- DFG/USBtoHDMI 轉換器
- 電源供應器
- 微型 HDMI 連接線
- 快速入門指南

www.theimagingsource.com



Headquarters: The Imaging Source Europe GmbH
Überseetor 18, 28217 Bremen, Germany
Phone: +49 421 33591-0

North & South America: The Imaging Source, LLC
6926 Shannon Willow Rd, S 400, Charlotte, NC 28226, USA
Phone: +1 704-370-0110

Asia Pacific: The Imaging Source Asia Co. Ltd.
2F, No. 8, Xinhua 1st Road,
Taipei City 114, Neihu District, Taiwan
Phone: +886 2-2792-3153

www.theimagingsource.com

All product and company names in this document may be trademarks and tradenames of their respective owners and are hereby acknowledged.

The Imaging Source Europe GmbH cannot and does not take any responsibility or liability for any information contained in this document. The source code presented in this document is exclusively used for didactic purposes. The Imaging Source does not assume any kind of warranty expressed or implied, resulting from the use of the content of this document or the source code.

The Imaging Source Company reserves the right to make changes in specifications, function or design at any time and without prior notice.

Last update: July 2020

Copyright © 2020 The Imaging Source Europe GmbH

All rights reserved. Reprint, also in parts, only allowed with permission of The Imaging Source Europe GmbH.

All weights and dimensions are approximate. Unless otherwise specified the lenses shown in the context of cameras are not shipped with these cameras.